ENVIRONMENTAL CHAMBER

CATALOGUE NO.: GEC-P40604


Environmental Chamber is used for conditioning of the samples before test. It is suitable for conditioning of all type of paper and board before and during testing.

Test Description:

A sample of a paper or board is conditioned and said to be in equilibrium with the conditioning atmosphere when it has been exposed to the atmosphere long enough for the result of two weighing of a test piece performed at and interval of at least 1 hour and not to differ by more than 0.25 per cent.

APPLICABLE STANDARDS:

TAPPI, SCAN P 2, ISO



Salient Features:

Complete with sealed unit and contact thermometer.

Digital display of Temperature and Humidity.

Global Engineering Corporation

Global Engineering Corporation



82,Ganga Sagar Scheme, 200 Feet Bypass, Sirsi Road, Jaipur-302021 Tel: +91-141-5179040,2470290, Fax: +91-141-5179040

Email: sales@gecpaper.com

Visit us at: www.gecpaper.com

© Copyright 2006 Global Engineering Corporation. Contact for the detail. Rev 00-0